

ENOCHS Treatment Tables

When budgets or room size are deciding factors in selecting examination room equipment or a full-featured table is not necessary, the ENOCHS #49 Jumbo and the #49 Economy Treatment Tables are perfect choices. Versatility in use and durability in construction, both tables are ideal for Diagnostic, Sports Medicine, Orthopedics or Physical Therapy.

Jumbo 49D Treatment Table

Dimensions: Height - 32"
Length - 72"
Width - 32"

One-piece Upholstered Top design is free of cross seams, decorative embossing and fabric folds that collect dirt and make cleaning difficult. Durable easy to clean Naugahyde® upholstery, covering 1½" of high-density foam, is manufactured to be wear, stain, mold, mildew, and oil resistant.

Paper Roll Holder is designed for easy loading of 18" or 21" wide paper roll.

A **32" Extra-wide Top** allows patient to feel secure during examination or treatment

Two 4" Deep Drawers with steel reinforced side guide suspension system to provide smooth discipline movement

Attractive **Formica™ base** for the rich look of solid wood

Steel Reinforced H-braced design providing strength and dimensional stability



ENOCHS is an ISO 9001:2000 and WOB Certified Company

Economy 48 Treatment Table

Dimensions: Height - 33"
Length - 72"
Width - 24"

One-piece Upholstered Top design is free of cross seams, decorative embossing and fabric folds that collect dirt and make cleaning difficult. Durable easy to clean Naugahyde® upholstery, covering 1½" of high-density foam, is manufactured to be wear, stain, mold, mildew, and oil resistant.

A **24" Wide Top** is ideal when space is limited














Fully Welded All-Steel Frame base is designed to provide strength and dimensional stability.

Paper Roll Holder is concealed and designed for easy loading of 18" or 21" wide paper roll



ENOCHS is an ISO 9001:2000 and WOB Certified Company

For additional information on the ENOCHS Treatment Tables and other fine examination room furniture call Customer Service at 800-428-2305, FAX to 317-580-2944, e-mail to enochs@enochsmed.com or write to P.O. Box 50559 Indianapolis, IN 46250

Feature	ENOCHS JUMBO 49/49D AND ECONOMY 48 TREATMENT TABLES SPECIFICATIONS				
<p>Color Selections</p> <p>Colors and textures represented here, are an approximate of true color under the limits of electronic viewing.</p> <p>Contact ENOCHS for actual color samples.</p> <p>Colors subject to change without notification</p>	<p>48 Base Colors</p>  <p>Pearl Grey</p>  <p>Cream White</p>  <p>Mocha</p>	<p>49 & 49D Base Colors</p>  <p>Grey</p>  <p>Beige</p>  <p>Natural Oak</p>	<p>Upholstery Top Colors for all Treatment Tables</p>  <p>Smoke</p>  <p>Windsor Wine</p>  <p>Cambridge Blue</p>  <p>Ninja</p>  <p>Earthen Brown</p>  <p>Fawn</p>  <p>Emerald</p>		
<p>Base Construction</p>	<p>JUMBO 49 and 49D Tables</p> <ul style="list-style-type: none"> • Steel reinforced H-brace design for strength and dimensional stability • Base finished in durable smooth surface Formica™ for fast efficient cleaning <p>ECONOMY 48 Table</p> <ul style="list-style-type: none"> • Fully welded all steel base designed for strength and durability • Powder coated metal surfaces for long service life • Metal surface has a durable smooth finish for quick efficient cleaning • Assembly required 				
<p>Top Construction</p>	<p>For All Treatment Tables</p> <ul style="list-style-type: none"> • One-piece seamless upholstered top design is free of cross-seams, decorative embossing and fabric folds for quick and effective cleaning • Covered with durable Naugahyde™ manufactured to be wear, stain, mold, mildew resistant • Padded with 1.5" high density foam for patient comfort • Upholstered sections securely attached to table frame, eliminating areas of contamination associated with removable top pads. • Top engineered for easy replacement if required 				
<p>Storage</p>	<p>Jumbo 49 Treatment Table No storage available</p> <p>Jumbo 49D Treatment Table Two 4" deep drawers total of 1.5 cubic feet of storage space</p> <ul style="list-style-type: none"> • High impact, rounded interior corners and edges for easy cleaning • Steel reinforced, side mounting drawer system providing smooth disciplined movement of drawers • Drawer front covered front, back and edges with durable easy to clean laminate <p>Economy 48 Table No storage available</p>				
<p>Specifications</p>	<p>JUMBO 49 and 49D Tables</p> <p>Length - 72" Width - 32" Height - 32"</p>		<p>ECONOMY 48 Table</p> <p>Length - 72" Width - 24" Height - 33"</p>		
<p>Warranty</p>	<p>4 – YEAR WARRANTY for all treatment tables</p>				
<p>Replacement Parts After Warranty Ends</p>	<p>ENOCHS strives to use standard commercial available parts in ENOCHS products. After warranty period ends, most parts can be purchased from ENOCHS or from local supply house which allows facility to secure the most competitive price for the parts.</p>				

ORDERING INFORMATION

ENOCHS can only be ordered through Authorized ENOCHS distributors. Call at 800-428-2305 or e-mail enochs@enochsmed.com for the distributor in your area.

QUICK-SHIP * Treatment Tables

Product Code

49D Grey Base and Drawer Fronts , Smoke Upholstered Top	49D-GRSM
48 Pearl Grey Base, Smoke Upholstered Top	48-PGSM
48 Pearl Grey Base, Windsor Wine Upholstered Top	48-PGWW
48 Pearl Grey Base, Cambridge Blue Upholstered Top	48-PGCB

Other Popular Treatment Table Color Combinations **

Jumbo 49D Treatment Table with Drawers

49D Grey Base and Drawer Fronts, Cambridge Blue Upholstered Top	49-DGRCB
49D Grey Base and Drawer Fronts , Windsor Wine Upholstered Top	49-DGRWW
49D Grey Base and Drawer Fronts , Ninja Upholstered Top	49-DGRNJ
49D Beige Base and Drawer Fronts , Fawn Upholstered Top	49-DBEFN
49D Natural Oak Base and Drawer Fronts , Fawn Upholstered Top	49-DNOFN

49 Jumbo Treatment Table

49 Grey Base, Smoke Upholstered Top	49-GRSM
49 Grey Base, Cambridge Blue Upholstered Top	49-GRCB
49 Grey Base, Windsor Wine Upholstered Top	49-GRWW
49 Grey Base, Ninja Upholstered Top	49-GRNJ
49 Beige Base, Fawn Upholstered Top	49-BEFN
49 Natural Oak Base, Fawn Upholstered Top	49-NOFN

48 Economy Treatment Table

48 Pearl Grey Base, Smoke Upholstered Top	48-GRSM
48 Pearl Grey Base, Cambridge Blue Upholstered Top	48-GRCB
48 Pearl Grey Base, Windsor Wine Upholstered Top	48-GRWW
48 Cream White Base, Fawn Upholstered Top	48-CWFFN
48 Mocha Base, Smoke Upholstered Top	48-MBSM

* Quick-ship products are those commonly kept in the ENOCHS inventory. Stock levels fluctuate based on orders. Call 800-428-2305 for availability.

** Non-inventoried product. When ordered, product will be scheduled into normal monthly production.