



ENOCHS ENCORE

compare.

	ENCORE 4250	ENCORE 4200
Warranty:	4-YEAR	2-YEAR
Base Construction:	Fully welded all steel base	
Power System:	Electro-Mechanical	
Power Adjustments:	Height: 24" - 40" Back: horizontal - 80°	
Foot Control:	"Hands free", standard on every table, SELV	
Hand Control:	"Easy to read", optional accessory, add at anytime, SELV	
Couch Top:	100% seamless vacuum formed style, all steel frame construction	
Perineal Recess:	Moulded, no seams, 10" W x 2" D at center point	
Pelvic Tilt:	7°, 11°	
Leg Extension:	Extends out 20" from end of seat section	
Irrigation Basin:	4" Deep stainless steel basin	
Electrical Outlets:	Hospital grade, UL Listed, dual electrical outlets with dedicated circuit breaker	Not available on this model
Auto Return/Stop	Not available on this model	
Instrument Warmer: with ON/OFF Switch, Circuit Breaker & Reset	High efficiency warmer utilizing footend drawer	Not available on this model
Master ON/OFF Switch	Shuts off power to entire table	
Stirrups:	Multi-positional, lateral and length adjustment	
Storage:	Side drawer: 15.5" W x 15" D x 6" H Footend drawer: 15.5" W x 15" D x 4" H Total Storage Space: 1.4 cubic feet	
Paper Roll Holder:	Accommodates 18" or 21" paper widths	
Accessories:	hand control, knee crutches, padded arm board, drawer locks, restraint straps, IV pole, urology basin with screen and drain, headend accessory rails, footend accessory rails, caster base	
Mains Power Entry:	IEC connection with internal circuit breaker, detachable cord	
Duty Cycle:	Table: 1 minute ON / 9 minutes OFF Receptacles: Continuous operation Instrument Warming Drawer: Continuous operation	Table: 1 minute ON / 9 minutes OFF Receptacles: Not applicable Instrument Warming Drawer: Not applicable
Ratings:	Table: 115V 60 Hz 10 amps Receptacles: 115V 60 Hz Total receptacle load: 3.9 amps	Table: 115V 60 Hz 10 amps Receptacles: Not applicable Total receptacle load: Not applicable
Table Dimensions:	Width: 27", Height: 24" - 40", Length: 52" (extends to 72")	
Table Weight:	416 lbs.	
Patient Load:	500 lbs., evenly distributed on main couch	
Certifications:	UL Listed	