



ENOCHS EXCEL

compare.

	EXCEL 356	EXCEL 355	EXCEL 350
Warranty:	4-YEAR	4-YEAR	2-YEAR
Base Construction:	Fully welded all steel base		
Back Lift Mechanism:	Horizontal to 80° by spring assisted counter balance		
Couch Top:	One-piece upholstered top, steel reinforced		
Perineal Recess:	Not available on this model		
Pelvic Tilt:	11°	Not available on this model	
Leg Extension:	Extends out 18" from end of seat section		
Irrigation Basin:	4" Deep basin		
Electrical Outlets:	Hospital grade, UL Listed, dual electrical outlets		Not available on this model
Instrument Warmer: with ON/OFF Switch, Circuit Breaker & Reset	High efficiency warmer utilizes top footend drawer	Not available on this model	
Stirrups:	Multi-positional, lateral and length adjustment		
Storage:	3 Side drawers: 21.5" W x 15" D x 4" H 2 Footend drawers: 15.5" W x 15" D x 4" H Total Storage Space: 3.32 cubic feet		
Foot Step:	Extra long surface engineered for 600 lb. patient load		
Paper Roll Holder:	Accommodates 18" or 21" paper widths		
Accessories:	knee crutches, drawer locks, stainless steel basin, restraint straps, IV pole, urology basin with drain and screen, footend accessory rails, caster base		
Mains Power Entry:	IEC Connection, internal circuit breaker with reset	IEC Connection, internal circuit breaker with reset	Not applicable
Duty Cycle:	Receptacles: Continuous operation Instrument Warming Drawer: Continuous operation	Receptacles: Continuous operation Instrument Warming Drawer: Not applicable	Receptacles: Not applicable Instrument Warming Drawer: Not applicable
Ratings:	Table: 115V 60 Hz 10 amps Receptacles: 115V 60Hz Total receptacle load: 9 amps		Table: Not applicable Receptacles: Not applicable Total Load: Not applicable
Table Dimensions:	Width: 27", Height: 34", Length: 55" (extends to 73")		
Table Weight:	268 lbs.		
Patient Load:	600 lbs., evenly distributed on main couch		
Certifications:	by UL to Standard 60601-1, CAN/CSA C22.3 No. 601-1		