

# Specifications for the ENOCHS Horizon 2250 High-Low Treatment Table

<b>Warranty:</b>	4-YEAR
<b>Base Construction:</b>	Fully welded all steel base
<b>Power System:</b>	Electro-Mechanical provides smooth and quiet movement
<b>Power Adjustments:</b>	<b>Height:</b> 23" - 39"
<b>Foot Control:</b>	"Hands free", standard on every table, SELV
<b>Hand Control:</b>	Not available on this model
<b>Couch Top:</b>	One-piece upholstered top
<b>Perineal Recess:</b>	Not available on this model
<b>Pelvic Tilt:</b>	Not available on this model
<b>Leg Extension:</b>	Not available on this model
<b>Irrigation Basin:</b>	Not available on this model
<b>Electrical Outlets: with Circuit Breaker</b>	Not available on this model
<b>Instrument Warmer: with ON/OFF Switch, Circuit Breaker &amp; Reset</b>	Not available on this model
<b>Stirrups:</b>	Not available on this model
<b>Storage:</b>	Not available on this model <b>Total Storage Space:</b> n/a
<b>Paper Roll Holder:</b>	Accommodates 18" or 21" paper widths
<b>Accessories:</b>	Not available for this model
<b>Mains Power Entry: with Circuit Breaker</b>	IEC connection, internal circuit breaker with reset
<b>Duty Cycle:</b>	<b>Table:</b> 1 minute ON / 9 minutes OFF <b>Receptacles:</b> n/a <b>Instrument Warming Drawer:</b> n/a
<b>Ratings:</b>	<b>Table:</b> 115V 60 Hz 10 amps <b>Receptacles:</b> n/a <b>Total Receptacle Load:</b> n/a
<b>Table Dimensions:</b>	<b>Width:</b> 32", <b>Height:</b> 23" - 39", <b>Length:</b> 72"
<b>Table Weight:</b>	350 lbs.
<b>Patient Load:</b>	500 lbs. evenly distributed on main couch
<b>Certifications:</b>	By UL to Standard 60601-1, CAN/CSA C22.3 No. 601-1