

ENOCHS RESOLUTE

compare.



	RESOLUTE 6250	RESOLUTE 6220	RESOLUTE 6210
Warranty:	4-YEAR		2-YEAR
Base Construction:	Fully welded all steel base		
Power System:	Electro-mechanical		
Programmable:	Yes	Programmable with optional hand control	Not available on this model
Memory Positions:	Four	Four with optional hand control	Not available on this model
Attainable Positions:	Chair, Contoured Chair, Trendelenburg, OB/GYN, Lithotomy, Urology, Flat		
Power Adjustments:	Height: 24" to 40" Back: horizontal to 80°		Tilt: horizontal to 30° Legrest: vertical to horizontal
Auto-Return:	Attainable as a programmed position	Attainable as a programmed position with optional hand control	Not available on this model
Foot Control:	Easy touch, "hands free" foot control Safely Extra-Low Voltage (SELV) and IPX1		
Hand Control:	Standard, Programmable, Safety Extra-Low Voltage (SELV)	Optional accessory that may be added in the field	Not available on this model
Hand Control Connection Port:	Connection ports are standard on both sides of table's moveable column Hand control moves with the table, eliminating over-stretched cord		Not available on this model
Upholstered Couch Top:	One-piece style, no cross seams, embossing, or folded fabric. Securely fastened to all steel frame - no Velcro™		
Headend Accessory Rails:	Stainless steel rails, mounted on both sides of the table		
Irrigation Basin:	2" deep stainless steel basin, integrated safety switch prohibits operation of foot section when basin is extended		
Safety Lockout:	Accessible from both sides of table, shuts off power to all articulating sections of table, but retains power to receptacles		Not available on this model
Electrical Receptacles:	Four hospital grade electrical receptacles (2 per side)		Not available on this model
Receptacle Overload Protection:	Dedicated circuit breaker with quick reset - no fuses		Not available on this model
Stirrups:	Positive lock multi-position stirrups, length and lateral adjustment, conceal when not in use		
Paper Roll Holder:	Accommodates 18" or 21" paper widths		
Accessories:	Knee crutches, chair arms, hand control (for 6220 only), padded arm board, restraint straps, IV pole, urology basin with screen and drain, footend accessory rails, caster base, universal clamp, multi-positional clamp		
Mains Power Entry:	Master ON/OFF, IEC connection, overload protection: integrated circuit breaker with quick reset - no fuses		
Duty Cycle:	1 minute ON / 9 minutes OFF		
Mode of Operation:	Table: Intermittent operation Receptacles: Continuous operation		Table: Intermittent operation Receptacles: Not applicable
Ratings:	Table: 115V 60Hz 10 amps Receptacles: 115V 60 Hz Total receptacle load: 3.9 amps		Table: 115V 60Hz 10 amps Receptacles: Not applicable Total receptacle load: Not applicable
Table Dimensions:	Length: 71", extends to 83" Width top: 30", seat tapers to 25" at headend Headrest: 19" wide, tapers to 15" Legrest: 19 wide, tapers to 18"		
Table Weight:	485 lbs.		
Patient Load:	500 lbs., evenly distributed on main couch		
Certifications:	UL Standard 60601-1 CAN/CSA C22.2 No. 601.1		