

Specifications for the ENOCHS REGENCY 456 OB/GYN Examination Table

Warranty:	4-YEAR
Base Construction:	Fully welded all steel base
Back Lift Mechanism:	Horizontal to 80° by spring assisted counter balance
Couch Top:	100% Seamless vacuum formed style All steel frame construction
Perineal Recess:	Moulded, no seams, 10" W x 2" D at center point
Pelvic Tilt:	7°, 11°
Leg Extension:	Extends out 20" from end of seat section
Irrigation Basin:	4" Deep basin
Electrical Outlets:	Dual electrical outlets, hospital grade, UL Listed
Instrument Warmer: with ON/OFF Switch, Circuit Breaker & Reset	High efficiency warmer utilizes top footend drawer
Stirrups:	Multi-positional, lateral and length adjustment
Storage:	3 Side Drawers: 21.5" W x 15" D x 4" H 2 Footend Drawers: 15.5" W x 15" D x 4" H Total Storage Space: 3.32 cubic feet
Foot Step:	Extra long surface engineered for 600 lb. patient load
Paper Roll Holder:	Accommodates 18" or 21" paper widths
Accessories:	knee crutches, arm board, drawer locks, stainless steel basin, restraint straps, IV pole, urology basin with drain and screen, headend accessory rails, footend accessory rails, caster base
Mains Power Entry: with Circuit Breaker	Standard feature
Duty Cycle:	Receptacles: Continuous Instrument Warming Drawer: Continuous
Ratings:	Table: 115V 60Hz 10 amps Receptacles: 115V 60Hz Total Receptacle Load: 9 amps
Table Dimensions:	Width: 27", Height: 34", Length: 52" (extends to 72")
Table Weight:	268 lbs.
Patient Load:	600 lbs. evenly distributed on main couch
Certifications:	By UL to Standard 60601-1, CAN/CSA C22.3 No. 601-1