

ENOCHS PED TABLE™

PEDIATRIC EXAMINATION TABLE

One-piece Upholstered Top design is free of cross seams, decorative embossing and fabric folds that collect dirt and make cleaning difficult. Durable easy to clean Naugahyde® upholstery, covering 1½" of high-density foam, is manufactured to be wear, stain, mold, mildew, and oil resistant.

Inch/centimeter Measure is mounted safely inside display case to eliminate cleaning issues.

Easy load **Paper Roll Holder** at head end of table accommodates 18" or 21" paper.

Table Base is a fully welded all steel construction for strength and long service life

Double Door Storage Cabinet with sliding half shelf provides for storage of bulk supplies



Dimensions

Overall Table
Height - 36"
Width - 50"
Depth - 24"

Ped Table Top
Length - 50"
Width - 18"
Back Guard - 6"
Front Guard - 4"

Side Cabinet
Height - 32"
Width - 25"
Depth - 18"

Display Case/Back Guard made of durable Tuffek® is manufactured to have high abrasion, chemical and impact resistance.

#51 Side Cabinet with two 4" deep drawers and double door storage area is an ideal base to hold the type of scale that works best for each individual facility.

Colorful characters and geometric shapes capture child's attention making visit fun and enjoyable.

Front Guard is 4" high and padded with ½" foam.

Total of seven 4" deep **Storage Drawers** with rounded internal corners and edges for quick easy cleaning

ENOCHS is an ISO 9001:2000 and WOB Certified Company



The Ped Table™ offers 8.56 cubic feet of storage space distributed among seven 4" deep drawers and a double door bulk storage area.

The drawers are constructed of high impact polystyrene with rounded corners and edges so cleaning is fast and efficient. Table is designed with a steel reinforced side guide mounting system for smooth non-binding drawer movement. Automatic stop eliminates accidental spillage of drawer contents.

Bulk storage area is engineered without a center bar which allows un-restricted access to stored items. A removeable sliding metal half-shelf allows the ability to store tall items along with other bulk supplies.
















Colorful characters and geometric shapes is what makes this pediatric table popular with the patient as well as the parents. The impact resistant plexiglass display case also acts as a secure back guard. The easy to read metal inch/centimeter measuring rule is mounting inside the display case to keep it clean and safely away from patient. The one-piece upholstered top covered in durable Naugahyde™ is designed for quick and efficient cleaning after each patient.



The Enoch's Ped Table™ is designed without an integrated scale so that each individual practice can choose the type of scale which best suits their specific needs. The separate scale concept eliminates downtime of the table if a scale malfunctions. A separate independent scale allows a facility to reduce overall costs by having the option of a central weighing station instead of a scale in every room.

The 51 Side Cabinet not only is an ideal unit to hold a scale but increases storage space with two 4" deep drawers and a double door bulk storage area with a sliding half shelf.

For additional information on the Enoch's Ped Table™ and other fine examination room furniture call Customer Service at 800-428-2305, FAX to 317-580-2944, e-mail to enochs@enochsmed.com or write to P.O. Box 50559 Indianapolis, IN 46250

Feature	ENOCHS PED TABLE AND SIDE CABINET																	
<p>Color Selections</p> <p>Colors and textures represented here, are an approximate of true color under the limits of electronic viewing.</p> <p>Contact ENOCHS for actual color samples.</p> <p>Colors subject to change without notification</p>	<p>Base Colors</p>  <p>Pearl Grey</p>  <p>Cream White</p>  <p>Mocha</p>	<p>Drawer/Door Front Colors</p>  <p>Grey</p>  <p>Beige</p>  <p>Natural Oak</p>	<p>Upholstery Top Colors</p>  <p>Smoke</p>  <p>Windsor Wine</p>  <p>Cambridge Blue</p>  <p>Ninja</p>  <p>Earthen Brown</p>  <p>Fawn</p>  <p>Emerald</p>															
<p>Base Construction of Ped 50 and 51 Cabinet</p>	<ul style="list-style-type: none"> • 18 gauge steel reinforced with heavier gauge steel for strength • Fully welded base and top construction for dimensional stability • Powder coated metal surfaces for durability and long service life • Base surfaces are smooth for easy efficient cleaning 																	
<p>Top Construction</p>	<p>Ped 50 Table Top</p> <ul style="list-style-type: none"> • One-piece seamless upholstered top design eliminates cross-stitching, decorative embossing and fabric folds for quick and effective cleaning • Padded with 1.5" high density foam for patient comfort • Covered in durable Naugahyde™ manufactured to be wear, sulfide stain, mold, mildew and oil resistant • Top securely attached to table, no removable parts • Convenient paper roll holder and cutter <p>Ped 51 Side Cabinet</p> <ul style="list-style-type: none"> • Seamless one-piece Formica® work surface 																	
<p>Display Case /Back Guard</p>	<ul style="list-style-type: none"> • Constructed of impact resistant Tuffek® plexiglass designed for abrasion and chemical resistance. • Encloses colorful characters and geometric shapes • Houses an easy to read inch/centimeter metal measuring rule • Functions as a unique 6" back guard. 																	
<p>Storage</p>	<p>Ped 50 Table</p> <p>Seven (7) 4" deep drawers offering a total of 5.2 cubic feet of drawer space Double door bulk storage compartment with 4.6 cubic feet of storage area</p> <p>Ped 51 Side Cabinet</p> <p>Two(2) 4" deep drawers offering a total of 1.5 cubic feet of drawer space Double door bulk storage compartment with 4.6 cubic feet of storage area</p> <p>Information Common to both Units</p> <ul style="list-style-type: none"> • Drawers are high impact and manufactured with round interior corners and edges for easy cleaning • Steel reinforced side guide mounting system for smooth disciplined movement of drawer • Drawer fronts are covered front, back and edges with durable easy to clean laminate • Bulk area has a removable sliding half-shelf to allow storage of tall items • Designed without center support bar allowing full unrestricted access into bulk storage area • Locks for drawers and doors available for additional fee 																	
<p>Specifications</p>	<table border="0"> <tr> <td>Table top:</td> <td>Overall table:</td> <td>Side Cabinet:</td> </tr> <tr> <td>Length - 50"</td> <td>Length - 50"</td> <td>Length - 25"</td> </tr> <tr> <td>Width - 18"</td> <td>Width - 24"</td> <td>Depth - 18"</td> </tr> <tr> <td>Back guard - 6"</td> <td>Height - 36"</td> <td>Height - 32"</td> </tr> <tr> <td>Front guard - 4"</td> <td></td> <td></td> </tr> </table>			Table top:	Overall table:	Side Cabinet:	Length - 50"	Length - 50"	Length - 25"	Width - 18"	Width - 24"	Depth - 18"	Back guard - 6"	Height - 36"	Height - 32"	Front guard - 4"		
Table top:	Overall table:	Side Cabinet:																
Length - 50"	Length - 50"	Length - 25"																
Width - 18"	Width - 24"	Depth - 18"																
Back guard - 6"	Height - 36"	Height - 32"																
Front guard - 4"																		
<p>Warranty</p>	<p>4 – YEAR WARRANTY</p>																	
<p>Replacement Parts After Warranty Ends</p>	<p>ENOCHS strives to use standard commercial available parts in ENOCHS products. After warranty period ends, most parts can be purchased from ENOCHS or from local supply house which allows facility to secure the most competitive price for the parts.</p>																	

ORDERING INFORMATION

ENOCHS can only be ordered through Authorized ENOCHS distributors. Call at 800-428-2305 or e-mail enochs@enochsmed.com for the distributor in your area.

Popular Pediatric Examination Table Color Combinations**	Product Code
Grey Drawer Fronts, Pearl Grey Base, Smoke Upholstered Top	50-GRPGSM
Grey Drawer Fronts, Pearl Grey Base, Windsor Wine Upholstered Top	50-GRPGWW
Grey Drawer Fronts, Pearl Grey Base, Cambridge Blue Upholstered Top	50-GRPGCB
Grey Drawer Fronts, Pearl Grey Base, Ninja Upholstered Top	50-GRPGNJ
Grey Drawer Fronts, Pearl Grey Base, Emerald Upholstered Top	50-GRPGEM
Beige Drawer Fronts, Cream White Base, Cambridge Blue Upholstered Top	50-BECWCB
Beige Drawer Fronts, Cream White Base, Windsor Wine Upholstered Top	50-BECWWW
Beige Drawer Fronts, Cream White Base, Fawn Upholstered Top	50-BECWFN
Natural Oak Drawer Fronts, Cream White Base, Fawn Upholstered Top	50-NOCWFN
Grey Drawer Fronts, Mocha Base, Smoke Upholstered Top	50-GRMBSM

Pediatric Side Cabinet Color Combinations **	Product Code
Grey Drawer Fronts, Pearl Grey Base	51-GRPG
Beige Drawer Fronts, Cream White Base	51-BECW
Natural Oak Drawer Fronts, Cream White Base	51-NOCW
Grey Drawer Fronts, Mocha Base	51-GRMB

** Non-inventoried product. When ordered, product will be scheduled into normal monthly production.

For additional information on the ENOCHS Ped Table™ and other fine examination room furniture call Customer Service at 800-428-2305, FAX to 317-580-2944, e-mail to enochs@enochsmed.com or write to P.O. Box 50559 Indianapolis, IN 46250