

ENOCHS POWER 2000-01

HIGH-LOW TREATMENT TABLE

Attainable Positions: Flat

Power System: Electro Hydraulic

Adjustments by "Easy Touch" Foot Control:

Height - 23" to 39"

Dimensions: Length - 72"

Width - 32"

Height - 23" to 39"

Accessories: Restraint straps

One-piece upholstered top design eliminates cross seams, decorative embossing and crevices that collect dirt and make cleaning difficult. Covered in durable, easy to clean Naugahyde® upholstery manufactured to be wear, stain, mold, mildew and oil resistant.

Paper Roll Holder is concealed at head end of table and accommodates 18" or 21" paper.

Paper Cutter located at headend of table.

"Easy Touch" Foot Control adjusts table's height with reliable hands-free action. Pedal is covered with non-slip tread and utilizes a natural side-to-side movement of operator's foot.

Padded with 1½" **High-density Foam** for patient comfort.

The 16" Range of Height adjustment allows patient to be positioned to a comfortable height for staff.

Table Rating: 120 VAC 60 Hz 12 A

Patient Load: 400# evenly distributed

Table is UL Listed

ENOCHS is an ISO 9001:2000 and WOB Certified Company



"Easy Touch" foot control allows smooth adjustment of table's height. Use of the foot control reduces incidents of cross contamination possible with hand-held controls. Operator utilizes the natural side-by-side movement of the foot on the control pedal to adjust the patient to desired height for the examination or treatment.





The extra wide 32" top of the Power 2000 makes it an ideal high-low table for general treatment, physical therapy and casting. The one-piece upholstered top is free of cross seams and decorative embossing which promotes quick efficient cleaning after each patient. Top is padded with 1½" of high-density foam for patient comfort.



Table quickly and smoothly adjusts in height from a low 23" to a full table height of 39". The table provides low access for disabled and geriatric patients in wheelchairs and those using walking aids. The low access point also reduces back stress or injury of staff assisting patient on or off the table. The height adjustment range of 16" allows physician or staff member to position patient at convenient working height.

For additional information on the ENOCHS Power 2000-01 and other fine examination room furniture call Customer Service at 800-428-2305, FAX to 317-580-2944, e-mail to enochs@enochsmed.com or write to P.O. Box 50559 Indianapolis, IN 46250

| Feature | ENOCHS POWER 2000-01 SPECIFICATIONS | |
|---|---|--|
| <p>Color Selections</p> <p>Colors and textures represented here, are an approximate of true color under the limits of electronic viewing.</p> <p>Contact ENOCHS for actual color samples.</p> <p>Colors subject to change without notification</p> | <p>Base Colors</p>  <p>Cream White Pearl Grey</p> <p>Mocha</p> | <p>Upholstery Top Colors</p>  <p>Smoke Windsor Wine Cambridge Blue</p> <p>Earthen Brown Fawn Ninja Emerald</p> |
| Base Construction | <ul style="list-style-type: none"> • 18 gauge steel reinforced with heavier gauge steel for strength • Fully welded base for dimensional stability • Powder coated metal surfaces for durability and long service life • Column style base with smooth surfaces for easy efficient cleaning | |
| Top Construction | <ul style="list-style-type: none"> • 100% one-piece upholstered top design is free of cross-seams, decorative embossing and fabric folds making cleaning quick and effective. • Padded with 1½" high-density foam for patient comfort • Covered in durable Naugahyde® manufactured to be wear, stain, mold, mildew and oil resistant • Upholstered section is securely attached to top frame • Top engineered for easy replacement if required | |
| Power System | <ul style="list-style-type: none"> • Electro-Hydraulic | |
| Power Adjustment | <ul style="list-style-type: none"> • Table Height : 23" to 39" | |
| Controls | <ul style="list-style-type: none"> • <i>Hands Free</i> foot control for fast accurate positioning of patient • Foot control eliminates cross contamination possible with hand controls and enhances infection control procedures. | |
| Specifications | <ul style="list-style-type: none"> • Length - 71" • Width - 32" • Height - 23" to 39" • Table Rating: 120 VAC 60 Hz 12 A • Entire Product UL Listed | |
| Warranty | 2 – YEAR WARRANTY | |
| Replacement Parts After Warranty Ends | <p>ENOCHS strives to use standard commercial available parts in ENOCHS products. After warranty period ends, most parts can be purchased from ENOCHS or from local supply house which allows facility to secure the most competitive price for the parts.</p> | |

ORDERING INFORMATION

ENOCHS can only be ordered through Authorized ENOCHS distributors.
 Call 800-428-2305 or e-mail enochs@enochsmed.com for the distributor in your area.

Popular Power 2000-01 Color Combinations**

Product Code

| | |
|---|-------------|
| Pearl Grey Base, Smoke Upholstered Top | 2000-01PGSM |
| Pearl Grey Base, Windsor Wine Upholstered Top | 2000-01PGWW |
| Pearl Grey Base, Cambridge Blue Upholstered Top | 2000-01PGCB |
| Pearl Grey Base, Ninja Upholstered Top | 2000-01PGNJ |
| Pearl Grey Base, Emerald Upholstered Top | 2000-01PGEM |
| Cream White Base, Earthen Brown Upholstered Top | 2000-01CWEB |
| Cream White Base, Fawn Upholstered Top | 2000-01CWFN |
| Mocha Base, Smoke Upholstered Top | 2000-01MBSM |
| Mocha Base, Earthen Brown Upholstered Top | 2000-01MBEB |
| Mocha Base, Fawn Upholstered Top | 2000-01MBFN |

** Non-inventoried product. When ordered, product will be scheduled into normal monthly production.