

# ENOCHS POWER 6000-01

## PROCEDURAL TABLE

**Padded Headrest** has a double arm linkage that allows versatile adjustment for patient comfort.

**Paper Roll Holder** is concealed and designed for easy loading of 18" or 21" wide paper roll.

**Head end Accessory Rails** on both sides of table allows quick easy attachment of accessories.

**One-piece upholstered top** design eliminates cross-seams, decorative embossing and crevices that collect dirt and make cleaning difficult. Easy to clean Naugahyde® upholstery, covering 1.5" of high density foam, is manufactured to be wear, stain, mold, mildew and oil resistant.

**Electrical Outlets** and the **Auto Return/Stop** controls are located on both sides of the table for convenient access.

**Master ON/OFF Power Switch** shuts off all power to foot control and auto return feature to prevent accidental movement of table during examination or procedure.

Hospital grade and **UL Listed Dual Electrical Outlets** – 2 per side.

**Auto Return Stop** overrides the Auto Return command and immediately stops table's decent.

**Auto Return** returns table to lowest position after examination or procedure.

**Table Rating:** 115VAC, 60 Hertz, 12 ampere  
**Patient Load:** 400# evenly distributed  
Table is UL Listed

**Attainable Positions:** Chair, Contoured Chair, Proctologic, Trendelenburg, OB/GYN, Lithotomy, Urology and Flat

**Power System:** Electro-Hydraulic

**Adjustments by "Easy Touch" Foot Control:**

Height – 26" to 42"  
Back – horizontal to 87° chair  
Tilt – horizontal to 42° Trendelenburg  
Leg – vertical to horizontal

**Dimensions:** Length - 71" Extended length - 80"  
Width - Top-27" at seat, tapers to 25" at headend  
Headrest - 19" wide, tapers to 15"  
Footrest - 19" wide, tapers to 17"  
Legrest - 19" wide

**Accessories:** Knee crutches, padded armboard, restraint straps, IV Pole, urology basin with screen and drain, chair arms, footend accessory rails and caster base



Positive lock **Multi-position Stirrups** with length and lateral adjustment, conceal when not in use.

2" deep **Stainless Steel Irrigation Basin** conveniently slides in and out from under seat section.

**Padded Leg rest** adjusts to support patient's legs during exam or procedure. An integrated safety feature will not allow leg section to operate if irrigation basin is not in full stored position under the seat section. Padded Legrest easily and quickly converts into a Kneerest/Footstep without awkward handles or levers.

**"Easy Touch" Foot Control** adjusts table's height, back, tilt and leg/foot rest with reliable *hands-free* action. Pedals are covered with non-slip tread and utilize a natural side-to-side movement of operator's foot.

## ENOCHS is an ISO 9001:2000 and WOB Certified Company



### OB/GYN Position

Multi-position stirrups adjust quickly and easily to desired length and into a standard or extra-wide angle. The system's positive locking design eliminates the dangerous bar overhang. Power 6000's ability to adjust in tilt as well as height makes it a perfect table for a variety of examination, treatment or procedures.



### Proctological Position

The 6000 seat section adjusts to various degrees from horizontal to 42° to provide desired tilt for examination, procedure or treatment. The fully padded legrest adjusts quickly to any one of four kneerest positions without awkward release handles or levers. The tables' height is adjustable from 26" to 42" to smoothly bring patient to a height convenient for the physician.














### Lithotomy & Urology Positions

Assisting in maintaining the abduction of thigh in Lithotomy Position, the use of the knee crutch accessory allows versatility of positioning. Other accessories for the 6000 include Urology basin with screen, drain and hose, IV pole, padded armboard, chair arms, restraint straps, foot end accessory rails and caster base.

For additional information on the Enoch's Power 6000-01 and other fine examination room furniture call Customer Service at 800-428-2305, FAX to 317-580-2944, e-mail to [enochs@enochsmed.com](mailto:enochs@enochsmed.com) or write to P.O. Box 50559 Indianapolis, IN 46250

## Feature

## ENOCHS 6000-01 PROCEDURAL TABLE SPECIFICATIONS

<p><b>Color Selections</b></p> <p>Colors and textures represented here, are an approximate of true color under the limits of electronic viewing.</p> <p>Contact ENOCHS for actual color samples.</p> <p>Colors subject to change without notification</p>	<p style="text-align: center;"><b>Base Colors</b></p> <div style="display: flex; justify-content: space-around;"> <div style="text-align: center;">               Cream White         </div> <div style="text-align: center;">               Pearl Grey         </div> </div> <div style="text-align: center; margin-top: 20px;">               Mocha         </div>	<p style="text-align: center;"><b>Upholstery Top Colors</b></p> <div style="display: flex; justify-content: space-around;"> <div style="text-align: center;">               Smoke         </div> <div style="text-align: center;">               Windsor Wine         </div> <div style="text-align: center;">               Cambridge Blue         </div> </div> <div style="display: flex; justify-content: space-around; margin-top: 20px;"> <div style="text-align: center;">               Earthen Brown         </div> <div style="text-align: center;">               Fawn         </div> <div style="text-align: center;">               Ninja         </div> <div style="text-align: center;">               Emerald         </div> </div>				
<p><b>Base Construction</b></p>	<ul style="list-style-type: none"> <li>• 18 gauge steel reinforced with heavier gauge steel for support and dimensional stability</li> <li>• Fully welded base and top construction</li> <li>• State-of-the-art powder coating on all metal surfaces for durability and long service life</li> </ul>					
<p><b>Top Construction</b></p>	<ul style="list-style-type: none"> <li>• Steel frame for strength, durability and long service life</li> <li>• One-piece upholstered top design eliminates cross-seams, decorative embossing and fabric folds</li> <li>• Design addresses important issues of cross contamination and infection control</li> <li>• Padded with 1.5" of high density foam</li> <li>• Covered in durable Naugahyde® manufactured to be wear, stain, mold, mildew and oil resistant</li> <li>• Top securely attaches to table frame -- no removable parts</li> </ul>					
<p><b>Attainable Position</b></p>	<p>Chair, Contoured Chair, OB/GYN, Lithotomy, Urology, Trendelenburg, Pediatric and Flat</p>					
<p><b>Control Mechanism</b></p>	<p><i>Hands Free "Easy Touch"</i> foot control to maintain best possible infection control procedures and reduce cross contamination</p>					
<p><b>Adjustment Features by Power Control</b></p>	<ul style="list-style-type: none"> <li>• Table Height : 26" to 42"</li> <li>• Back : Horizontal to 87°</li> <li>• Tilt: Horizontal to 42°</li> <li>• Legrest: Vertical to Horizontal</li> </ul>					
<p><b>One-piece Legrest Pad</b></p>		<p>The 6000-01 features a one-piece padded leg section offering a large work surface for leg and foot procedures while providing patient with comfortable support. The unit is padded with 1.5" of high density foam, covered in durable Naugahyde to match main couch and is securely attached to the foot section frame.</p>				
<p><b>Power System</b></p>	<p>Electro-Hydraulic</p>					
<p><b>Electrical</b></p>	<p>UL listed hospital grade electrical outlets located on both sides of table for easy access. Total of 4 outlets, 2 per side.</p>					
<p><b>Headrest</b></p>	<p>Articulating headrest</p>					
<p><b>Stirrups</b></p>	<ul style="list-style-type: none"> <li>• Multi-position, three lateral positions</li> <li>• Adjustment of shoe distance from end of table</li> <li>• Shoe and bar move simultaneously for length adjustment, eliminating dangerous bar overhang</li> <li>• Vinyl coated shoe for comfort and ease of cleaning</li> </ul>					
<p><b>Accessory Rails</b></p>	<ul style="list-style-type: none"> <li>• Headend accessory rails are a standard feature</li> <li>• footend accessory rails optional</li> </ul>					
<p><b>Accessories</b></p>	<table style="width: 100%; border: none;"> <tr> <td style="vertical-align: top; width: 50%;"> <ul style="list-style-type: none"> <li>• universal clamp (part# 6109)</li> <li>• tri-clamp (part# 6110)</li> <li>• footend accessory rails (part# 6112)</li> <li>• knee crutches<sup>†</sup> (part# 6122)</li> <li>• padded armboard* (part# 6124)</li> </ul> </td> <td style="vertical-align: top; width: 50%;"> <ul style="list-style-type: none"> <li>• restraint straps (part# 6126)</li> <li>• IV Pole* (part# 6130)</li> <li>• urology pan with drain, screen &amp; hose (part# 6137)</li> <li>• chair arms (part# 6140)</li> <li>• caster base (part# 6150)</li> </ul> </td> </tr> <tr> <td style="text-align: center; font-size: small;">* requires #6109 or 6110 clamp</td> <td style="text-align: center; font-size: small;">† requires # 6112 footend accessory rails</td> </tr> </table>		<ul style="list-style-type: none"> <li>• universal clamp (part# 6109)</li> <li>• tri-clamp (part# 6110)</li> <li>• footend accessory rails (part# 6112)</li> <li>• knee crutches<sup>†</sup> (part# 6122)</li> <li>• padded armboard* (part# 6124)</li> </ul>	<ul style="list-style-type: none"> <li>• restraint straps (part# 6126)</li> <li>• IV Pole* (part# 6130)</li> <li>• urology pan with drain, screen &amp; hose (part# 6137)</li> <li>• chair arms (part# 6140)</li> <li>• caster base (part# 6150)</li> </ul>	* requires #6109 or 6110 clamp	† requires # 6112 footend accessory rails
<ul style="list-style-type: none"> <li>• universal clamp (part# 6109)</li> <li>• tri-clamp (part# 6110)</li> <li>• footend accessory rails (part# 6112)</li> <li>• knee crutches<sup>†</sup> (part# 6122)</li> <li>• padded armboard* (part# 6124)</li> </ul>	<ul style="list-style-type: none"> <li>• restraint straps (part# 6126)</li> <li>• IV Pole* (part# 6130)</li> <li>• urology pan with drain, screen &amp; hose (part# 6137)</li> <li>• chair arms (part# 6140)</li> <li>• caster base (part# 6150)</li> </ul>					
* requires #6109 or 6110 clamp	† requires # 6112 footend accessory rails					
<p><b>On-Off Switch</b></p>	<p>Master ON/OFF Power Switch shuts off all power to foot control and auto return feature to prevent accidental movement of table during examination or procedure.</p>					

Feature	<b>ENOCHS 6000 PROCEDURAL TABLE SPECIFICATIONS</b>	
<b>Auto Return and Stop</b>	Auto Return returns table to lowest position after examination or procedure. The Auto Return Stop overrides the Auto Return command and immediately stops table's decent. Start and stop controls are accessible from either side of table	
<b>Specifications</b>	<ul style="list-style-type: none"> <li>• Length - 71" extends to 80"</li> <li>• Maximum Width - 27"</li> <li>• 2 -1/2" deep stainless steel irrigation basin</li> </ul>	<ul style="list-style-type: none"> <li>• Table rating: 115 VAC 60 cycle 12 amps</li> <li>• UL Listed hospital grade receptacles</li> <li>• Entire table is UL Listed</li> </ul>
<b>Warranty</b>	<b>4 – YEAR WARRANTY</b>	
<b>Replacement Parts After Warranty Ends</b>	Standard commercial available parts are used in ENOCHS products. After warranty period ends, parts can be purchased from ENOCHS or from local supply house which allows facility to secure the most competitive price for the parts.	

## ORDERING INFORMATION

ENOCHS can only be ordered through Authorized ENOCHS distributors. Call at 800-428-2305 or e-mail [enochs@enochsmed.com](mailto:enochs@enochsmed.com) for the distributor who serves your area.

### **Power 6000-01 Procedural Tables Popular Color Combinations\*\***

### **Product Code**

Pearl Grey Base, Smoke Upholstered Top	6000-01PGSM
Pearl Grey Base, Windsor Wine Upholstered Top	6000-01PGWW
Pearl Grey Base, Cambridge Blue Upholstered Top	6000-01PGCB
Pearl Grey Base, Ninja Upholstered Top	6000-01PGNJ
Cream White Base, Fawn Upholstered Top	6000-01CWFN
Pearl Grey Base, Emerald Upholstered Top	6000-01PGEM
Cream White Base, Cambridge Blue Upholstered Top	6000-01CWCB
Cream White Base, Windsor Wine Upholstered Top	6000-01CWWW
Cream White Base, Ninja Upholstered Top	6000-01CWNJ
Cream White Base, Earthen Brown Upholstered Top	6000-01CWEB
Cream White Base, Emerald Upholstered Top	6000-01CWEM
Mocha Base, Smoke Upholstered Top	6000-01MBSM
Mocha Base, Windsor Wine Upholstered Top	6000-01MBWW
Mocha Base, Earthen Brown Upholstered Top	6000-01MBEB
Mocha Base, Fawn Upholstered Top	6000-01MBFN

\*\* Non-inventoried product. When ordered, product will be scheduled into normal monthly production.

For additional information on the ENOCHS Power 6000-01 and other fine examination room furniture call Customer Service at 800-428-2305, FAX to 317-580-2944, e-mail to [enochs@enochsmed.com](mailto:enochs@enochsmed.com) or write to P.O. Box 50559 Indianapolis, IN 46250