

# Enochs Stools, Blood Drawing Chair & Waste Receptacles

An examination room or lab is completed with ENOCHS stools, waste receptacles, side and blood drawing chairs. In the tradition of exceptional quality, these units are among the strongest and best performing on the market.

ENOCHS stools have fully padded seats covered in durable easy to care for Naugahyde®. Stools are available with manual or pneumatic bases and are engineered to assure smooth stable movement.

All waste receptacles feature ENOCHS EXTENDED LIFE CONSTRUCTION® offering strength, durability and top performance of a state-of-the-art powder coated finish. Designs range from a jumbo capacity unit, suitable for high traffic facilities with high volume of disposable paper products to smaller units perfect for moderately used examination rooms.

## Color Selections:

See ENOCHS Color Collection for available colors



### 106 Pneumatic Stool with Back

offers style and performance. The five leg, polished aluminum base provides smooth stable movement. Stools height adjusts quickly and easily by a hand lever located under the seat. Padded back features a 3" vertical adjustment.

**Adjustment: 20" to 30"**

**101 Pneumatic Stool** is designed with a five leg base. The comfortable padded seat is upholstered in durable Naugahyde®.

The 101 adjusts smoothly and easily by control lever under the seat. The 2" carpet casters allow quiet and controlled movement.

**Adjustment: 20" to 30"**



**96 Operator's Stool with Back** has a threaded stem adjustment. The five leg, powder coated, tubular steel base is fully welded and features a sturdy foot ring. Padded back features 3" vertical adjustment. Stool is covered in durable, easy to clean Naugahyde®. The 2" ball bearing casters provide smooth quiet movement.

**Adjustment: 18.5" to 25"**



**92 Operator's Stool** is designed with a fifth leg to increase stability. State-of-the-art powder coated tubular steel base is fully welded and features an attractive foot ring. The padded seat is upholstered in durable easy to maintain Naugahyde® and adjusts in height by seat rotation.

**Adjustment: 18.5" to 25"**



**201 Blood Drawing Chair** incorporates an extra-wide concept in its design that provides non-restricted access for the patient.

Padded arm boards have full lateral and a 25" to 38" height adjustment.

The stainless steel side cabinet features a seamless work surface and one drawer for efficient storage. Cabinet is adjustable in height to meet technician's needs and is easily changed for left or right side access.

**45" W x 21" D x 30" H**

**200 Side Chair** is part of the 201 and features fully welded chrome base with a polypropylene seat and back. Designed for strength, style and ease of cleaning, this unit is perfect in any room.

**21" W x 21" D x 30" H**



**112 Jumbo Step-on Waste Receptacle** has a 10.5 gallon capacity providing efficient handling of the large amount of disposable paper goods generated in high traffic medical facilities. The single "step-on" action lifts the lid for fast, hands-free disposal of trash. A steel sanitary guard surrounds the removable internal plastic container, assuring waste containment.

This hefty unit has an all steel, fully welded base with a state-of-the-art powder coated finish and accented by laminate panel and lid.

**Capacity: 10.5 gallon**  
**17" W x 17" D x 30" H**

**109 Waste Receptacle** has a fully welded steel base finished in durable state-of-the-art powder coating and accented with a laminate panel to match ENOCHS cabinets or tables. Standard with an open concave removable stainless steel lid as pictured on the 108 receptacle, the 109 is pictured here with an optional stainless steel flip top lid.

**Capacity: 7 gallon**  
**10" W x 10" D x 24" H**



**108 Waste Receptacle** is a fully welded steel base unit finished in durable powder coating. The open concave stainless steel lid holds liner bag in place and lifts off quickly for waste removal. Unit may be ordered with an optional stainless steel flip top lid as pictured on the 109 receptacle.

**Capacity: 7 gallon**  
**10" W x 10" D x 24" H**

# ENOCHS

EXAMINING ROOM FURNITURE

**800-428-2305**

## THE ENOCHS STOOLS, WASTE RECEPTACLES, BLOOD DRAWING AND SIDE CHAIR

### GENERAL OVERVIEW

ENOCHS offers various stools, waste receptacles and a blood drawing chair and side chair to complete any examination or lab area in a private practice, clinic or hospital.

ENOCHS stools have fully padded seats covered in durable easy to care for Naugahyde®. Stools are available with manual or pneumatic height adjustment and are engineered with five-leg bases to assure smooth stable movement. Stools are also available with padded backs.

All waste receptacles feature ENOCHS EXTENDED LIFE CONSTRUCTION® offering strength, durability and top performance of a state-of-the-art powder coated finish. Designs range from a jumbo capacity unit, suitable for high traffic facilities with high volume of disposable paper products to smaller units perfect for moderately used examination rooms.

ENOCHS Blood Drawing Chair features an extra wide design that provides non-restricted access for the patient. Two adjustable arm boards and a side cabinet complete the unit.

### SPECIFICATIONS

OPERATOR STOOLS	# 92	# 96	# 101	# 106
<b>Base type:</b>	tubular steel	tubular steel	cast star base	cast star base
<b>Number of legs:</b>	5 leg	5 leg	5 leg	5 leg
<b>Foot ring:</b>	yes	yes	no	no
<b>Base finish:</b>	power coated	power coated	polished aluminum	polished aluminum
<b>Adjustment type:</b>	threaded stem	threaded stem	pneumatic	pneumatic
<b>Adjustment height :</b>	18.5" to 25"	18.5" to 25"	20" to 30"	20" to 30"
<b>Back support:</b>	no	yes	no	yes
<b>Back adjustment:</b>	N/A	3" vertical	N/A	3" vertical
<b>Caster type :</b>	2" ball bearing, hard surface caster	2" ball bearing, hard surface caster	2" carpet caster	2" carpet caster
<b>Seat dimensions:</b>	14" W x 2.5" H	14" W x 2.5" H	14" W x 2.5" H	14" W x 2.5" H
<b>Warranty :</b>	4 year	4 year	4 year	4 year
<b>Shipping wt:</b>	43*	43*	47*	47*
<b>Actual wt:</b>	16	21	13	18
<b>Shipping box for base:</b>	5 cubic feet	5 cubic feet	5 cubic feet	5 cubic feet
<b>Shipping box for seat:</b>	included with base	1 cubic feet	included with base	1 cubic feet
<b>Where stools are used:</b>	Specialties: All		Placement: Exam or treatment rooms	

#### Color selection:

Colors and textures represented here are an approximate of the true color under the limits of electronic viewing. Contact ENOCHS for actual color samples.

Colors subject to change without notification



Cambridge Blue



Windsor Wine



Ninja



Emerald



Desert Sand



Smoke



Earthen Brown



Fawn



English Brown



Black

WASTE RECEPTICALS	# 108	# 109	# 112
<b>Base type:</b>	Fully welded steel construction	Fully welded steel construction	Fully welded steel construction
<b>Decorative accent:</b>	N/A	laminate front panel	laminate front panel and lid
<b>Base finish:</b>	power coated , treated for resistance to rust	power coated, treated for resistance to rust	power coated, treated for resistance to rust
<b>Trash access:</b>	Concave stainless top with non-restrictive access	Concave stainless top with non-restrictive access	Step on action lifts lid for disposal of waste
<b>Options:</b>	Stainless steel flip top lid	Stainless steel flip top lid	N/A
<b>Special features:</b>	N/A	N/A	Steel sanitary guard surrounding inner plastic receptacle to assure trash containment
<b>Capacity:</b>	7 gallon	7 gallon	10.5 gallon
<b>Dimensions:</b>	10" W x 10" D x 24" H	10" W x 10" D x 24" H	17" W x 17" D x 30" H
<b>Warranty :</b>	4 year	4 year	4 year
<b>Shipping wt:</b>	20*	20*	68
<b>Actual wt:</b>	13	13	60
<b>Shipping carton:</b>	2 cubic feet	2 cubic feet	7 cubic feet
<b>Where receptacles are used:</b>	Specialties: All Placement: Exam or treatment rooms		

**Color selection:**

Colors and textures represented here are an approximate of the true color under the limits of electronic viewing. Contact ENOCHS for actual color samples.

Colors subject to change without notification



<b>BLOOD DRAWING &amp; SIDE CHAIR</b>	<b>#201 BLOOD DRAWING CHAIR</b>	<b># 200 SIDE CHAIR</b>
<b>Base type:</b>	Fully welded steel construction	Fully welded steel construction
<b>Base finish:</b>	Chrome	Chrome
<b>Seat &amp; back finish :</b>	Smooth surface, polypropylene	Smooth surface, polypropylene
<b>Armboards:</b>	Two padded armboards, covered in durable Naughahyde®	N/A
<b>Armboard adjustment:</b>	lateral: full height: 25" to 38"	N/A
<b>Side cabinet:</b>	Steel cabinet with one drawer, designed to mount under armboard on left or right-hand side depending on room configuration	N/A
<b>Dimension:</b>	45" W x 21" D x 30" H	21" W x 21" D x 30" H
<b>Warranty :</b>	4 year	4 year
<b>Shipping wt:</b>	96*	96*
<b>Actual wt:</b>	56	18
<b>Shipping carton:</b>	11 cubic feet	11 cubic feet
<b>Where receptacles are used:</b>	Specialties: All Placement: lab area for blood draws	Specialties: All Placement: Areas for side chairs
<b>Color selection:</b>	Black	Black

\* Dimensional weight is required to calculate shipping charges