

ENOCHS RESOLUTE 6250

Programmable Procedural Table

PATIENT LOAD: 500 LBS.
evenly distributed on the table's main couch



TABLE FEATURES:

- 4-YEAR WARRANTY
- Programmable Hand Control with 4 memory positions
- "Easy Touch" foot control for hands free table positioning
- Power Adjustment of Back, Tilt, Foot and Height
- Safety Lockout to prevent accidental movement
- Circuit Breakers for overload protection-no fuses
- Electrical Receptacles, Safety Lockouts and Hand Control Ports located on table's movable column for convenience and easy access
- Cords plugged into table's receptacles, move up and down with the table, eliminating over-stretched cords
- Master ON/OFF
- Multi-Position Positive Lock Stirrups
- Stainless Steel Irrigation Basin
- Stainless Steel Headend Accessory Rails



Positioning includes chair, contoured chair, Trendelenburg, Lithomy, OBGYN and flat. A safety anti-pinch system stops movement of table if a pinch situation exists between headrest and floor.



Power adjustment of height (24"-40"), back (flat to 80°), tilt (flat to 30°) and foot sections (vertical to horizontal) allows the 6250 to be positioned for numerous exams and procedures.