Specifications for the ENOCHS ENCORE 4200 High-Low Examination Table

Warranty:	2-YEAR
Base Construction:	Fully welded all steel base
Power System:	Electro-Mechanical provides smooth and quiet movement
Power Adjustments:	Height: 24"- 40" Back: horizontal - 80 °
Foot Control:	"Hands free", standard on every table, SELV
Hand Control:	"Easy to read", optional accessory, add at anytime, SELV
Couch Top:	100% Seamless vacuum formed style, all steel frame construction
Perineal Recess:	Moulded, no seams, 10" W x 2" D at center point
Pelvic Tilt:	7°, 11°
Leg Extension:	Extends out 20" from end of seat section
Irrigation Basin:	4" Deep stainless steel basin
Electrical Outlets: with Circuit Breaker	Not available on this model
Instrument Warmer: with ON/OFF Switch, Circuit Breaker & Reset	Not available on this model
Stirrups:	Multi-positional, lateral and length adjustment
Storage:	Side Drawer: 15.5" W x 15" D x 6" H Footend Drawer: 15.5" W x 15" D x 4" H Total Storage Space: 1.4 cubic feet
Paper Roll Holder:	Accommodates 18" or 21" paper widths
Accessories:	hand control, knee crutches, padded arm board, drawer locks, restraint straps, IV pole, urology basin with screen and drain, headend accessory rails, footend accessory rails, caster base
Mains Power Entry: with Circuit Breaker	IEC connection, internal circuit breaker with reset
Duty Cycle:	Table: 1 minute ON / 9 minutes OFFReceptacles: n/aInstrument Warming Drawer: n/a
Ratings:	Table : 115V 60 Hz 10 amps Receptacles : n/a Total Receptacle Load : n/a
Table Dimensions:	Width: 27", Height: 24" - 40", Length: 52" (extends to 72")
Table Weight:	416 lbs.
Patient Load:	500 lbs. evenly distributed on main couch
Certifications:	By UL to Standard 60601-1, CAN/CSA C22.3 No. 601-1



800-428-2305

www.enochsmed.com enochs@enochsmed.com