

ENOCHS ENCORE 4510

BASIC LOW-ACCESS EXAMINATION TABLE

PATIENT LOAD: 700 LBS.
evenly distributed on the table's main couch



TABLE FEATURES:

- 2-YEAR Warranty
- Patient Load: 700 lbs.
- Low access: 18.5" Maximum height: 38"
- Power Back adjusts from horizontal to 80°
- Easy to clean 100% Seamless Vacuum Formed Top with additional padding on edges for patient comfort
- 11° Pelvic Tilt
- Stainless Steel Footend Drawer
- Stainless Steel Irrigation Basin operates independent of Leg Extension
- Multi-position Positive Lock Stirrups with length and lateral adjustment
- Overload Protection: Circuit Breakers - NO FUSES
- Fully Welded All Steel Base Construction - NO PLASTIC
- Classified by UL to ANSI/AAMI ES60601-1 (2005) + AMD 1 (2012), CAN/CSA-C22.2 No. 60601-1 (2014)



Multi-position stirrups, a padded pull-out leg extension and a separate pull-out Stainless Steel irrigation basin are all features of the 4510 designed to reduce time and increase efficiency in exams or procedures.



The 2" deep recess, moulded into the seat section, provides more comfortable positioning during perineal exams or procedures. The 11° pelvic tilt with release handles located on both sides of the table is a standard feature.



A Grab Rail System accessory for ENCORE tables provides assistance to patients when getting on or off the table. The unique design of ENOCHS' all steel grab rails allows staff to quickly remove one or both units when not needed and easily reinstall when required.