## **ENOCHS ENCORE** 4520

LOW-ACCESS EXAMINATION TABLE

**PATIENT LOAD: 700 LBS.** 

evenly distributed on the table's main couch



## **TABLE FEATURES:**

- 4-YEAR Warranty
- · Patient Load: 700 lbs.
- · Low access: 18.5" Maximum height: 38"
- Power Back adjusts from horizontal to 80°
- Easy to clean 100% Seamless Vacuum Formed Top with additional padding on edges for patient comfort
- 11° Pelvic Tilt
- Field Reversible Side Drawer
- Stainless Steel Footend Drawer
- Stainless Steel Irrigation Basin operates independent of Leg Extension
- Multi-position Positive Lock Stirrups with length and lateral adjustment
- Overload Protection: Circuit Breakers NO FUSES
- Fully Welded All Steel Base Construction NO PLASTIC
- Classified by UL to ANSI/AAMI ES60601-1 (2005) + AMD 1 (2012), CAN/CSA-C22.2 No. 60601-1 (2014)



Facilities create their perfect table by selecting ADD-ON OPTIONS they want, not forced to buy features they do not need. Options of hospital grade dual electrical receptacles with dedicated circuit breaker, instrument warmer with easy access ON/OFF switch and a standard non-programmable SELV hand control are all available to order with a table.



The ENOCHS 4520 is ADA compliant offering patients easy, no barrier access to table. The power back adjusts to any angle from horizontal to 80° and allows couch surface to be flat or in a chair position for patient transfers. A Grab Rail System to provide assistance for the patient is available as an accessory. The "Easy Touch" SELV foot control comes with the 4520 allowing hands free adjustment of table. A non-programmable hand control may be chosen as an ADD-ON OPTION or purchased at a later date.