ENOCHS ENCORE 4550

PROGRAMMABLE LOW-ACCESS EXAMINATION TABLE

PATIENT LOAD: 700 LBS.

evenly distributed on the table's main couch



TABLE FEATURES:

- 4-YEAR Warranty
- Patient Load: 700 lbs.
- Programmable Hand Control with 3 memory positions
- Low access: 18.5" Maximum height: 38"
- Power Back adjusts from horizontal to 80°
- Easy to clean 100% Seamless Vacuum Formed Top with additional padding on edges for patient comfort
- 11° Pelvic Tilt
- Field Reversible Side Drawer
- Stainless Steel Footend Drawer
- Stainless Steel Irrigation Basin operates independent of Leg Extension
- Multi-position Positive Lock Stirrups with length and lateral adjustment
- Overload Protection: Circuit Breakers NO FUSES
- Fully Welded All Steel Base Construction NO PLASTIC
- Classified by UL to ANSI/AAMI ES60601-1 (2005) + AMD 1 (2012), CAN/CSA-C22.2 No. 60601-1 (2014)



Working end of the table features multiposition stirrups and a padded pullout leg extension designed to move independent of the stainless steel irrigation basin for increased efficiency. The stainless steel footend drawer provides storage and accommodates the instrument warmer when ordered as an ADD-ON OPTION.



The ENOCHS 4550 is ADA compliant with safe, easy, no barrier access onto table's 18.5" low seat height. The power back adjusts to any angle from horizontal to 80° allowing couch surface to be flat or in a chair position for confident patient transfer. A Grab Rail System for patient assistance is available as an accessory.