

Specifications for the ENOCHS Horizon 2250 High-Low Treatment Table

Warranty:	4-YEAR
Base Construction:	Fully welded all steel base
Power System:	Electro-Mechanical provides smooth and quiet movement
Power Adjustments:	Height: 23" to 39"
Foot Control:	"Hands-free", standard on every table, SELV
Hand Control:	Not available on this model
Couch Top:	One-piece upholstered top
Perineal Recess:	Not available on this model
Pelvic Tilt:	Not available on this model
Leg Extension:	Not available on this model
Irrigation Basin:	Not available on this model
Electrical Outlets: with Circuit Breaker	Not available on this model
Instrument Warmer: with ON/OFF Switch, Circuit Breaker & Reset	Not available on this model
Stirrups:	Not available on this model
Storage:	Not available on this model Total Storage Space: N/A
Paper Roll Holder:	Accommodates 18" or 21" paper widths
Accessories:	Not available for this model
Mains Power Entry: with Circuit Breaker	IEC connection, internal circuit breaker with reset
Duty Cycle:	Table: 1 minute ON / 9 minutes OFF Receptacles: N/A Total receptacle load: N/A
Ratings:	Table: 115V 60 Hz 10 amps Receptacles: N/A Total receptacle load: N/A
Table Dimensions:	Width: 32", Height: 23" - 39", Length: 72"
Table Weight:	350 lbs.
Patient Load:	500 lbs. evenly distributed on main couch
Certifications:	Classified by UL to UL60601-1, IEC/EN 60601-1, CAN/CSA-C22.2 No. 601.1