ENOCHS REGENCY



compare.

	REGENCY 456	REGENCY 450
Warranty:	4-YEAR	2-YEAR
Base Construction:	Fully welded all steel base	
Back Lift Mechanism:	Horizontal to 80° by spring assisted counter balance	
Couch Top:	100% Seamless vacuum formed style, all steel frame construction	
Perineal Recess:	Moulded, no seams, 10" W x 2" D at center point	
Pelvic Tilt:	7°, 11°	
Leg Extension:	Extends out 20" from end of seat section	
Irrigation Basin:	4" Deep basin	
Electrical Outlets:	Hospital grade, UL Listed, dual electrical outlets	Not available on this model
Instrument Warmer: with ON/OFF Switch, Circuit Breaker & Reset	High efficiency warmer utilizes top footend drawer	Not available on this model
Stirrups:	Multi-positional, lateral and length adjustment	
Storage:	3 Side drawers: $21.5''$ W x $15''$ D x $4''$ H 2 Footend drawers: $15.5''$ W x $15''$ D x $4''$ H Total Storage Space: 3.32 cubic feet	
Foot Step:	Extra long surface engineered for 600 lb. patient load	
Paper Roll Holder:	Accommodates 18" or 21" paper widths	
Accessories:	knee crutches, padded arm board, drawer locks, stainless steel basin, restraint straps, IV pole, urology basin with screen and drain, headend accessory rails, footend accessory rails, caster base	
Mains Power Entry:	IEC connection, internal circuit breaker with reset	Not applicable
Duty Cycle:	Receptacles: Continuous operation Instrument Warming Drawer: Continuous operation	Receptacles: Not applicable Instrument Warming Drawer: Not applicable
Ratings:	Table: 115V 60 Hz 10 amps Receptacles: 115V 60 Hz Total receptacle load: 9 amps	Table: Not applicable Receptacles: Not applicable Total receptacle load: Not applicable
Table Dimensions:	Width: 27", Height: 34", Length 52" (extends to 72")	
Table Weight:	268 lbs.	
Patient Load:	600 lbs. evenly distributed on main couch	
Certifications:	Classified by UL to UL60601-1, IEC/EN 60601-1, CAN/CSA-C22.2 No. 601.1	