ENOCHS RESOLUTE 6220

General Procedural Table

PATIENT LOAD: 500 LBS.

evenly distributed on the table's main couch



TABLE FEATURES:

- 4-YEAR WARRANTY
- "Easy Touch" Foot Control for hands free adjustment of table
- Electrical Receptacles, 2 per side. Overload Protection: quick reset Circuit Breaker — no fuses
- Safety Lockout, accessible from either side of table, prevents accidental movement of table during exam or procedure, but maintains power to electrical receptacles
- Hand Control Connection Port, standard on both sides of table, allows hand control accessory to be added at any time
- Receptacles, Safety Lockout and Hand Control Ports located on table's moveable column for convenience and easy access
- Items plugged into side of table travel with the height of the table, eliminating over-stretched cords
- Adjustable Padded Headrest
- Multi-Position Positive Lock Stirrups
- 2" deep Stainless Steel Irrigation Basin
- Stainless Steel Headend Accessory Rails



Power adjustment of height (24"–40"), back (flat – 80°), tilt (flat – 30°) and foot sections (vertical – horizontal) allows table to be positioned for various exams and procedures including chair, contoured chair, Trendelenburg, Lithomy, OBGYN and flat.



Quick access to electrical receptacles and the safety lockout switch on the side of the table saves time. The innovative Lockout System disables the hand and foot control to eliminate accidental table movement during exams or procedures but retains power to table's receptacles.